

# FCS Express<sup>7</sup> Cytometry




## Transform Your Data Into Results

FCS Express<sup>TM</sup> **7** provides data transformation tools in an easy-to-use interface, enabling any scientist access to powerful tools for understanding high dimensional data sets. Drag and drop any transformation on data grids, heat maps, and plots to apply!

### *t*SNE

- Map high dimensional data in 2D plots
- Identify populations with colors in density plots or heat maps



### *K* Means


- Perform cluster analysis
- Plot newly created parameters against each other or existing ones



### *Integration*

Run your scripts directly in FCS Express to:

- Add new parameters
- Create new transformations



### *Parameter Math*

Use mathematical formulas to:

- Create new parameters
- Change existing parameters

### *Principal Component Analysis*

Use this orthogonal transformation to create new variables, known as principal components

Reshape your analysis today. Get pricing and your FREE 30-day trial at:  
[www.denovosoftware.com](http://www.denovosoftware.com)

